

For Research Use Only

AIF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67791-3-PBS



Basic Information

Catalog Number:

67791-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12400

GenBank Accession Number:

BC111065

GeneID (NCBI):

9131

UNIPROT ID:

O95831

Full Name:

apoptosis-inducing factor,
mitochondrion-associated, 1

Calculated MW:

609 aa, 66 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D9D1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50508-1, AIF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67791-2-PBS. Detection antibody: 67791-3-PBS. Standard: Ag12400. Range: 3.125-100 ng/mL.