

For Research Use Only

PKMYT1 Monoclonal antibody, PBS Only



Catalog Number: 67806-1-PBS

Featured Product

Basic Information

Catalog Number:

67806-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30849

GenBank Accession Number:

NM_004203

GeneID (NCBI):

9088

UNIPROT ID:

Q99640

Full Name:

protein kinase, membrane associated tyrosine/threonine 1

Calculated MW:

55 kDa

Observed MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

1D12D3

Applications

Tested Applications:

WB, Indirect ELISA, IHC

Species Specificity:

Human, mouse, rat

Background Information

Protein kinase, membrane associated tyrosine/threonine (PKMYT1) is a member of the Wee family of protein kinases. PKMYT1 acts as a negative regulator of the cell cycle by inactivating the CDK1-cyclinB complex through phosphorylation of Tyr14/Tyr15 and preventing cell cycle entering into mitosis at the G2/M transition. In mammalian somatic cells, PKMYT1 affects the assembly of the Golgi apparatus and endoplasmic reticulum during telophase of mitosis by inhibiting cyclin activities. Abnormal expression of PKMYT1 is closely associated with multiple types of cancers. PKMYT1 protein was mainly localized in the cytoplasm. (PMID: 31695486, PMID: 32801781)

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67806-1-Ig (PKMYT1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67806-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67806-1-Ig (PKMYT1 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67806-1-PBS in a different storage buffer formulation.



WB result of PKMYT1 antibody (67806-1-Ig; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PKMYT1 transfected MCF-7 cells. This data was developed using the same antibody clone with 67806-1-PBS in a different storage buffer formulation.