

For Research Use Only

PKMYT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67806-4-PBS



Basic Information

Catalog Number: 67806-4-PBS	GenBank Accession Number: NM_004203	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9088	CloneNo.: 1G11F7
Source: Mouse	UNIPROT ID: Q99640	
Isotype: IgG2b	Full Name: protein kinase, membrane associated tyrosine/threonine 1	
Immunogen Catalog Number: AG30849	Calculated MW: 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50951-2, PKMYT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67806-4-PBS. Detection antibody: 67806-5-PBS. Standard:Ag30849. Range: 0.781-100 ng/mL.