

For Research Use Only

SMCR7L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67808-4-PBS



Basic Information

Catalog Number:

67808-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13775

GenBank Accession Number:

BC002587

GeneID (NCBI):

54471

UNIPROT ID:

Q9NQG6

Full Name:

Smith-Magenis syndrome
chromosome region, candidate 7-like

Calculated MW:

463 aa, 51 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E1C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50714-2, SMCR7L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67808-2-PBS. Detection antibody: 67808-4-PBS. Standard: Ag13775. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50714-3, SMCR7L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67808-5-PBS. Detection antibody: 67808-4-PBS. Standard: Ag13775. Range: 0.391-100 ng/mL.