

For Research Use Only

SMCR7L Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67808-5-PBS



Basic Information

Catalog Number: 67808-5-PBS	GenBank Accession Number: BC002587	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 54471	CloneNo.: 3A12E10
Source: Mouse	UNIPROT ID: Q9NQG6	
Isotype: IgG1	Full Name: Smith-Magenis syndrome chromosome region, candidate 7-like	
Immunogen Catalog Number: AG13775	Calculated MW: 463 aa, 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50714-3, SMCR7L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67808-5-PBS. Detection antibody: 67808-4-PBS. Standard: Ag13775. Range: 0.391-100 ng/mL.