

For Research Use Only

DAPK1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67815-2-PBS

Basic Information

Catalog Number: 67815-2-PBS	GenBank Accession Number: BC113660	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1612	CloneNo.: 1B10F6
Source: Mouse	UNIPROT ID: P53355	
Isotype: IgG1	Full Name: death-associated protein kinase 1	
Immunogen Catalog Number: AG29838	Calculated MW: 1430 aa, 160 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50195-1, DAPK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67815-2-PBS. Detection antibody: 67815-3-PBS. Standard: Ag29838. Range: 1.563-100 ng/mL.