

For Research Use Only

# DAPK1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67815-3-PBS

## Basic Information

Catalog Number:

67815-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29838

GenBank Accession Number:

BC113660

GeneID (NCBI):

1612

UNIPROT ID:

P53355

Full Name:

death-associated protein kinase 1

Calculated MW:

1430 aa, 160 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H7G3

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50195-1, DAPK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67815-2-PBS. Detection antibody: 67815-3-PBS. Standard: Ag29838. Range: 1.563-100 ng/mL.