

For Research Use Only

# HGS Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67818-2-PBS

## Basic Information

<b>Catalog Number:</b> 67818-2-PBS	<b>GenBank Accession Number:</b> BC003565	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 9146	<b>CloneNo.:</b> 2E10E5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O14964	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> hepatocyte growth factor-regulated tyrosine kinase substrate	
<b>Immunogen Catalog Number:</b> AG28610	<b>Calculated MW:</b> 86 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50245-1, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-2-PBS. Detection antibody: 67818-3-PBS. Standard: Ag28610. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50245-2, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-2-PBS. Detection antibody: 67818-5-PBS. Standard: Ag28610. Range: 1.563-100 ng/mL.