

For Research Use Only

# HGS Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67818-3-PBS

## Basic Information

Catalog Number:

67818-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28610

GenBank Accession Number:

BC003565

GeneID (NCBI):

9146

UNIPROT ID:

O14964

Full Name:

hepatocyte growth factor-regulated  
tyrosine kinase substrate

Calculated MW:

86 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2C9F6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50245-1, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-2-PBS. Detection antibody: 67818-3-PBS. Standard: Ag28610. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50245-3, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-3-PBS. Standard: Ag28610. Range: 3.125-100 ng/mL.