

For Research Use Only

HGS Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67818-5-PBS

Basic Information

Catalog Number: 67818-5-PBS	GenBank Accession Number: BC003565	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9146	CloneNo.: 1D5E6
Source: Mouse	UNIPROT ID: O14964	
Isotype: IgG2a	Full Name: hepatocyte growth factor-regulated tyrosine kinase substrate	
Immunogen Catalog Number: AG28610	Calculated MW: 86 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50245-2, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-2-PBS. Detection antibody: 67818-5-PBS. Standard:Ag28610. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50245-4, HGS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67818-4-PBS. Detection antibody: 67818-5-PBS. Standard:Ag28610. Range: 1.563-100 ng/mL.