

For Research Use Only

SLC25A42 Monoclonal antibody, PBS Only



Catalog Number: 67819-1-PBS

Basic Information

Catalog Number: 67819-1-PBS	GenBank Accession Number: BC045598	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 284439	CloneNo.: 2F7G10
Source: Mouse	UNIPROT ID: Q86VD7	
Isotype: IgG2b	Full Name: solute carrier family 25, member 42	
Immunogen Catalog Number: AG26379	Calculated MW: 318 aa, 35 kDa	
	Observed MW: 33-35 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



K-562 cells were subjected to SDS PAGE followed by western blot with 67819-1-Ig (SLC25A42 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67819-1-PBS in a different storage buffer formulation.