

For Research Use Only

C9orf72 Monoclonal antibody, PBS Only



Catalog Number: 67824-1-PBS

Basic Information

Catalog Number: 67824-1-PBS	GenBank Accession Number: BC020851	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 203228	CloneNo.: 2F9A2
Source: Mouse	UNIPROT ID: Q96LT7	
Isotype: IgG2a	Full Name: chromosome 9 open reading frame 72	
Immunogen Catalog Number: AG21080	Calculated MW: 481 aa, 54 kDa	
	Observed MW: 50 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, FC

Species Specificity:
Human, mouse, rat, pig

Background Information

C9ORF72 has a domain with polymorphic hexanucleotide repeat (GGGGCC). The C9ORF72-hexanucleotide repeat expansions have been recently identified as genetic markers in amyotrophic lateral sclerosis (ALS) and frontotemporal lobar degeneration (FTLD). The C9ORF72 repeat expansions may indicate a worse prognosis in ALS. Human C9ORF72 has some isoforms with MW 54-60 kDa and 25-30 kDa. Mouse C9orf72 has some isoforms with MW 50-60 kDa and 35 kDa.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67824-1-Ig (C9orf72 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67824-1-PBS in a different storage buffer formulation.



1×10^6 HEK-293T cells were intracellularly stained with 0.4 ug Anti-Human C9orf72 (67824-1-Ig, Clone:2F9A2) and CoraLite® 488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (red), or 0.4 ug Mouse IgG2a Isotype Control (C1.18.4) (65208-1-Ig, Clone: C1.18.4) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 67824-1-PBS in a



Neuro-2a cells were subjected to SDS PAGE followed by western blot with 67824-1-Ig (C9orf72 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67824-1-PBS in a different storage buffer formulation.