

For Research Use Only

SH3BP5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67841-5-PBS



Basic Information

Catalog Number: 67841-5-PBS	GenBank Accession Number: BC010123	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9467	CloneNo.: 3D8G8
Source: Mouse	UNIPROT ID: O60239	
Isotype: IgG1	Full Name: SH3-domain binding protein 5 (BTK-associated)	
Immunogen Catalog Number: AG30252	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50694-2, SH3BP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67841-4-PBS. Detection antibody: 67841-5-PBS. Standard:Ag30252. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50694-3, SH3BP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67841-2-PBS. Detection antibody: 67841-5-PBS. Standard:Ag30252. Range: 0.391-50 ng/mL.