

For Research Use Only

VAPA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67845-2-PBS



Basic Information

Catalog Number: 67845-2-PBS	GenBank Accession Number: BC002992	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 9218	CloneNo.: 1D2D8
Source: Mouse	UNIPROT ID: Q9POL0	
Isotype: IgG1	Full Name: VAMP (vesicle-associated membrane protein)-associated protein A, 33kDa	
Immunogen Catalog Number: AG7819	Calculated MW: 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50773-1, VAPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67845-2-PBS. Detection antibody: 67845-3-PBS. Standard:Ag7819. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP50773-2, VAPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67845-2-PBS. Detection antibody: 67845-4-PBS. Standard:Ag7819. Range: 1.563-200 ng/mL.