

For Research Use Only

EMCN Monoclonal antibody, PBS Only



Catalog Number: 67854-1-PBS

Basic Information

Catalog Number: 67854-1-PBS	GenBank Accession Number: BC017781	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51705	CloneNo.: 2D8G3
Source: Mouse	UNIPROT ID: Q9ULC0	
Isotype: IgG1	Full Name: endomucin	
Immunogen Catalog Number: AG31146	Calculated MW: 19 kDa, 27 kDa	
	Observed MW: 120 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HUVEC cells were subjected to SDS PAGE followed by western blot with 67854-1-Ig (EMCN antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67854-1-PBS in a different storage buffer formulation.