

For Research Use Only

GSTM1 Monoclonal antibody, PBS Only



Catalog Number: 67862-1-PBS

Basic Information

Catalog Number: 67862-1-PBS	GenBank Accession Number: BC024005	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2944	CloneNo.: 3C5C2
Source: Mouse	UNIPROT ID: P09488	
Isotype: IgG2b	Full Name: glutathione S-transferase mu 1	
Immunogen Catalog Number: AG3113	Calculated MW: 181 aa, 21 kDa	
	Observed MW: 26-28 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, chicken, mouse, rat, pig

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat liver tissue were subjected to SDS PAGE followed by western blot with 67862-1-Ig (GSTM1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67862-1-PBS in a different storage buffer formulation.