

For Research Use Only

# HSP47/SERPINH1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67863-2-PBS



## Basic Information

<b>Catalog Number:</b> 67863-2-PBS	<b>GenBank Accession Number:</b> BC014623	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 871	<b>CloneNo.:</b> 2F7E5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P50454	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)	
<b>Immunogen Catalog Number:</b> AG30999	<b>Calculated MW:</b> 46 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50781-1, HSP47/SERPINH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67863-2-PBS. Detection antibody: 67863-1-PBS. Standard: Ag30999. Range: 0.781-100 ng/mL.