

For Research Use Only

# Synaptophysin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67864-5-PBS



## Basic Information

Catalog Number:

67864-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG11803

GenBank Accession Number:

BC064550

GeneID (NCBI):

6855

UNIPROT ID:

P08247

Full Name:

synaptophysin

Calculated MW:

313 aa, 34 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A11B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50359-2, Synaptophysin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67864-4-PBS. Detection antibody: 67864-5-PBS. Standard: Ag11803. Range: 3.125-50 ng/mL.