

For Research Use Only

NEDD4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67865-3-PBS



Basic Information

Catalog Number: 67865-3-PBS	GenBank Accession Number: BC152452	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4734	CloneNo.: 2H3A6
Source: Mouse	UNIPROT ID: P46934	
Isotype: IgG1	Full Name: neural precursor cell expressed, developmentally down-regulated 4	
Immunogen Catalog Number: AG17124	Calculated MW: 1319 aa, 149 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50851-1, NEDD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67865-2-PBS. Detection antibody: 67865-3-PBS. Standard: Ag17124. Range: 0.391-25 ng/mL