

For Research Use Only

# NEDD4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67865-5-PBS



## Basic Information

Catalog Number:

67865-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG17124

GenBank Accession Number:

BC152452

GeneID (NCBI):

4734

UNIPROT ID:

P46934

Full Name:

neural precursor cell expressed,  
developmentally down-regulated 4

Calculated MW:

1319 aa, 149 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B10H5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50851-2, NEDD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67865-4-PBS. Detection antibody: 67865-5-PBS. Standard:Ag17124. Range: 6.25-100 ng/mL



Cytometric bead array standard curve of MP50851-3, NEDD4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67865-6-PBS. Detection antibody: 67865-5-PBS. Standard:Ag17124. Range: 6.25-100 ng/mL