

For Research Use Only

# TMED10 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67876-2-PBS



## Basic Information

Catalog Number:

67876-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7516

GenBank Accession Number:

BC001825

GeneID (NCBI):

10972

UNIPROT ID:

P49755

Full Name:

transmembrane emp24-like trafficking protein 10 (yeast)

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F11C2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50860-1, TMED10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67876-2-PBS. Detection antibody: 67876-3-PBS. Standard: Ag7516. Range: 6.25-100 ng/mL