

For Research Use Only

TMED10 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67876-4-PBS



Basic Information

Catalog Number:

67876-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7516

GenBank Accession Number:

BC001825

GeneID (NCBI):

10972

UNIPROT ID:

P49755

Full Name:

transmembrane emp24-like trafficking protein 10 (yeast)

Calculated MW:

25 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E6E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50860-2, TMED10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67876-4-PBS. Detection antibody: 67876-3-PBS. Standard: Ag7516. Range: 3.125-100 ng/mL.