

For Research Use Only

FMO4 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67881-2-PBS



Basic Information

Catalog Number: 67881-2-PBS	GenBank Accession Number: BC002780	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2329	CloneNo.: 1B5H4
Source: Mouse	UNIPROT ID: P31512	
Isotype: IgG1	Full Name: flavin containing monooxygenase 4	
Immunogen Catalog Number: AG7935	Calculated MW: 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50902-1, FMO4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67881-2-PBS. Detection antibody: 67881-3-PBS. Standard:Ag7935. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50902-2, FMO4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67881-4-PBS. Detection antibody: 67881-2-PBS. Standard:Ag7935. Range: 0.098-100 ng/mL.