

For Research Use Only

SLC25A36 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 67896-2-PBS

Basic Information

Catalog Number:

67896-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23894

GenBank Accession Number:

BC014064

GeneID (NCBI):

55186

UNIPROT ID:

Q96CQ1

Full Name:

solute carrier family 25, member 36

Purification Method:

Protein G Magarose purification

CloneNo.:

2A4F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50339-1, SLC25A36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67896-2-PBS. Detection antibody: 67896-3-PBS. Standard: Ag23894. Range: 0.098-12.5 ng/mL



Cytometric bead array standard curve of MP50339-3, SLC25A36 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67896-6-PBS. Detection antibody: 67896-2-PBS. Standard: Ag23894. Range: 0.098-100 ng/mL