

For Research Use Only

# OSBPL5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67913-3-PBS



## Basic Information

|                                             |                                                       |                                                                |
|---------------------------------------------|-------------------------------------------------------|----------------------------------------------------------------|
| <b>Catalog Number:</b><br>67913-3-PBS       | <b>GenBank Accession Number:</b><br>BC032646          | <b>Purification Method:</b><br>Protein A Magarose purification |
| <b>Size:</b><br>1 mg/ml                     | <b>GeneID (NCBI):</b><br>114879                       | <b>CloneNo.:</b><br>1B11E3                                     |
| <b>Source:</b><br>Mouse                     | <b>UNIPROT ID:</b><br>Q9H0X9                          |                                                                |
| <b>Isotype:</b><br>IgG3                     | <b>Full Name:</b><br>oxysterol binding protein-like 5 |                                                                |
| <b>Immunogen Catalog Number:</b><br>AG11554 | <b>Calculated MW:</b><br>811 aa, 91 kDa               |                                                                |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50959-1, OSBP.L5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67913-2-PBS. Detection antibody: 67913-3-PBS. Standard:Ag11554. Range: 0.781-100 ng/mL.