

For Research Use Only

TRAPA, SSR1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67925-3-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 67925-3-PBS | GenBank Accession Number: BC007710 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 6745 | CloneNo.: 1E4A4 |
| Source: Mouse | UNIPROT ID: P43307 | |
| Isotype: IgG2b | Full Name: signal sequence receptor, alpha | |
| Immunogen Catalog Number: AG30832 | Calculated MW: 34 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50817-1, TRAPA, SSR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67925-2-PBS. Detection antibody: 67925-3-PBS. Standard: Ag30832. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50817-2, TRAPA, SSR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67925-4-PBS. Detection antibody: 67925-3-PBS. Standard: Ag30832. Range: 3.125-100 ng/mL.