

For Research Use Only

TRAPA, SSR1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67925-6-PBS



Basic Information

Catalog Number: 67925-6-PBS	GenBank Accession Number: BC007710	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 6745	CloneNo.: 2A10G11
Source: Mouse	UNIPROT ID: P43307	
Isotype: IgG2a	Full Name: signal sequence receptor, alpha	
Immunogen Catalog Number: AG30832	Calculated MW: 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50817-4, TRAPA, SSR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67925-4-PBS. Detection antibody: 67925-6-PBS. Standard: Ag30832. Range: 3.125-100 ng/mL.