

For Research Use Only

FAM114A1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:67926-3-PBS



Basic Information

Catalog Number:

67926-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG17289

GenBank Accession Number:

BC040452

GeneID (NCBI):

92689

UNIPROT ID:

Q8IWE2

Full Name:

family with sequence similarity 114, member A1

Calculated MW:

563 aa, 61 kDa

Purification Method:

Protein G purification

CloneNo.:

4E10F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51192-2, FAM114A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67926-3-PBS. Detection antibody: 67926-4-PBS. Standard:Ag17289. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51192-3, FAM114A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67926-3-PBS. Detection antibody: 67926-2-PBS. Standard:Ag17289. Range: 0.781-100 ng/mL