

For Research Use Only

PSME1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67928-3-PBS



Basic Information

Catalog Number: 67928-3-PBS	GenBank Accession Number: BC007503	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5720	CloneNo.: 1B4E2
Source: Mouse	UNIPROT ID: Q06323	
Isotype: IgG1	Full Name: proteasome (prosome, macropain) activator subunit 1 (PA28 alpha)	
Immunogen Catalog Number: AG30549	Calculated MW: 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51008-1, PSME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67928-2-PBS. Detection antibody: 67928-3-PBS. Standard: Ag30549. Range: 0.781-100 ng/mL.