

For Research Use Only

# RON, MST1R Monoclonal antibody, PBS Only



Catalog Number: 67954-1-PBS

## Basic Information

Catalog Number:

67954-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28289

GenBank Accession Number:

BQ936677 NM\_002447

GeneID (NCBI):

4486

UNIPROT ID:

Q04912

Full Name:

macrophage stimulating 1 receptor (c-met-related tyrosine kinase)

Calculated MW:

152 kDa

Observed MW:

150-185 kDa

Purification Method:

Protein G purification

CloneNo.:

2D8D1

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



BxPC-3 cells were subjected to SDS PAGE followed by western blot with 67954-1-Ig (RON, MST1R antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67954-1-PBS in a different storage buffer formulation.