

For Research Use Only

PSMB5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67959-3-PBS

Basic Information

Catalog Number: 67959-3-PBS	GenBank Accession Number: BC057840	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5693	CloneNo.: 1B10G5
Source: Mouse	UNIPROT ID: P28074	
Isotype: IgG1	Full Name: proteasome (prosome, macropain) subunit, beta type, 5	
Immunogen Catalog Number: AG5828	Calculated MW: 263 aa, 28 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50285-1, PSMB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67959-2-PBS. Detection antibody: 67959-3-PBS. Standard: Ag5828. Range: 0.781-200 ng/mL.