

For Research Use Only

PSMB5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67959-5-PBS

Basic Information

Catalog Number:

67959-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5828

GenBank Accession Number:

BC057840

GeneID (NCBI):

5693

UNIPROT ID:

P28074

Full Name:

proteasome (prosome, macropain)

subunit, beta type, 5

Calculated MW:

263 aa, 28 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A3E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50285-2, PSMB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67959-2-PBS. Detection antibody: 67959-5-PBS. Standard:Ag5828. Range: 0.781-200 ng/mL



Cytometric bead array standard curve of MP50285-3, PSMB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67959-4-PBS. Detection antibody: 67959-5-PBS. Standard:Ag5828. Range: 0.781-200 ng/mL