

For Research Use Only

HADHB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67967-2-PBS



Basic Information

Catalog Number:

67967-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30315

GenBank Accession Number:

BC017564

GeneID (NCBI):

3032

UNIPROT ID:

P55084

Full Name:

hydroxyacyl-Coenzyme A
dehydrogenase/3-ketoacyl-
Coenzyme A thiolase/enoyl-
Coenzyme A hydratase (trifunctional
protein), beta subunit

Calculated MW:

51 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D7D4

Applications

Tested Applications:

IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50546-1, HADHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67967-1-PBS. Detection antibody: 67967-2-PBS. Standard:Ag30315. Range: 3.125-100 ng/mL.