

For Research Use Only

# HADHB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67967-4-PBS



## Basic Information

Catalog Number:

67967-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG30315

GenBank Accession Number:

BC017564

GeneID (NCBI):

3032

UNIPROT ID:

P55084

Full Name:

hydroxyacyl-Coenzyme A  
dehydrogenase/3-ketoacyl-  
Coenzyme A thiolase/enoyl-  
Coenzyme A hydratase (trifunctional  
protein), beta subunit

Calculated MW:

51 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B10F4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50546-2, HADHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67967-3-PBS. Detection antibody: 67967-4-PBS. Standard: Ag30315. Range: 6.25-100 ng/mL