

For Research Use Only

BCAS2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68003-3-PBS

Basic Information

Catalog Number: 68003-3-PBS	GenBank Accession Number: BC005285	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10286	CloneNo.: 1C1E2
Source: Mouse	UNIPROT ID: O75934	
Isotype: IgG1	Full Name: breast carcinoma amplified sequence	
Immunogen Catalog Number: AG21309	2	
	Calculated MW: 26 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50169-1, BCAS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68003-2-PBS. Detection antibody: 68003-3-PBS. Standard: Ag21309. Range: 3.125-100 ng/mL.