

For Research Use Only

SFXN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68024-2-PBS



Basic Information

Catalog Number: 68024-2-PBS	GenBank Accession Number: BC020517	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 94081	CloneNo.: 2C8D6
Source: Mouse	UNIPROT ID: Q9H9B4	
Isotype: IgG1	Full Name: sideroflexin 1	
Immunogen Catalog Number: AG2985	Calculated MW: 35.6 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50630-1, SFXN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68024-2-PBS. Detection antibody: 68024-3-PBS. Standard: Ag2985. Range: 0.098-50 ng/mL