

For Research Use Only

FADS2-Specific Monoclonal antibody, PBS Only (Capture)

Catalog Number:68026-2-PBS



Basic Information

Catalog Number: 68026-2-PBS	GenBank Accession Number: BC009011	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9415	CloneNo.: 1B10G4
Source: Mouse	UNIPROT ID: O95864	
Isotype: IgG2a	Full Name: fatty acid desaturase 2	
Immunogen Catalog Number: AG27715	Calculated MW: 445 aa, 49 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50824-1, FADS2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68026-2-PBS. Detection antibody: 68026-3-PBS. Standard: Ag27715. Range: 0.195-25 ng/mL