

For Research Use Only

Endophilin 1 Monoclonal antibody, PBS Only



Catalog Number: 68029-1-PBS

Basic Information

Catalog Number:

68029-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG24793

GenBank Accession Number:

BC032825

GeneID (NCBI):

6456

UNIPROT ID:

Q99962

Full Name:

SH3-domain GRB2-like 2

Calculated MW:

39.8 kDa

Observed MW:

38-40 kDa

Purification Method:

Protein A purification

CloneNo.:

1G4H1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat, pig, rabbit

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68029-1-Ig (SH3GL2 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68029-1-PBS in a different storage buffer formulation.