

For Research Use Only

# mitoNEET,CISD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:68030-2-PBS



## Basic Information

Catalog Number: 68030-2-PBS	GenBank Accession Number: BC007043	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55847	CloneNo.: 1C9B6
Source: Mouse	UNIPROT ID: Q9NZ45	
Isotype: IgG2b	Full Name: CDGSH iron sulfur domain 1	
Immunogen Catalog Number: AG8560	Calculated MW: 108 aa, 12 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50573-1, mitoNEET, CISD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68030-2-PBS. Detection antibody: 68030-3-PBS. Standard: Ag8560. Range: 0.781-100 ng/mL