

For Research Use Only

RPLP0 Monoclonal antibody, PBS Only



Catalog Number: 68056-1-PBS

Basic Information

Catalog Number:

68056-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG1829

GenBank Accession Number:

BC009867

GeneID (NCBI):

6175

UNIPROT ID:

P05388

Full Name:

ribosomal protein, large, PO

Calculated MW:

34 kDa

Observed MW:

34 kDa

Purification Method:

Protein G purification

CloneNo.:

2A2E1

Applications

Tested Applications:

WB, Indirect ELISA, IHC

Species Specificity:

Human, Mouse, Rat

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68056-1-Ig (RPLP0 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68056-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 68056-1-Ig (RPLP0 antibody) at dilution of 1:400 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68056-1-PBS in a different storage buffer formulation.