

For Research Use Only

FATP2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68074-2-PBS



Basic Information

Catalog Number:

68074-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5390

GenBank Accession Number:

BC057770

GeneID (NCBI):

11001

UNIPROT ID:

O14975

Full Name:

solute carrier family 27 (fatty acid transporter), member 2

Calculated MW:

567 aa, 65 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G7G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50832-1, SLC27A2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68074-1-PBS. Detection antibody: 68074-2-PBS. Standard:Ag5390. Range: 12.5-200 ng/mL