

For Research Use Only

RAB12 Monoclonal antibody, PBS Only



Catalog Number: 68079-1-PBS

Basic Information

Catalog Number: 68079-1-PBS	GenBank Accession Number: NM_001025300	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 201475	CloneNo.: 2G11C2
Source: Mouse	UNIPROT ID: Q6IQ22	
Isotype: IgG2a	Full Name: RAB12, member RAS oncogene family	
Immunogen Catalog Number: AG31627	Calculated MW: 27KD	
	Observed MW: 27-30 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68079-1-Ig (RAB12 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68079-1-PBS in a different storage buffer formulation.