

For Research Use Only

C1QBP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68084-2-PBS



Basic Information

Catalog Number: 68084-2-PBS	GenBank Accession Number: BC013731	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 708	CloneNo.: 1A8B2
Source: Mouse	UNIPROT ID: Q07021	
Isotype: IgG2a	Full Name: complement component 1, q subcomponent binding protein	
Immunogen Catalog Number: AG19773	Calculated MW: 282 aa, 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50719-1, C1QBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68084-2-PBS. Detection antibody: 68084-3-PBS. Standard:Ag19773. Range: 0.391-100 ng/mL.