

For Research Use Only

# C1QBP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68084-3-PBS



## Basic Information

Catalog Number:

68084-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19773

GenBank Accession Number:

BC013731

GeneID (NCBI):

708

UNIPROT ID:

Q07021

Full Name:

complement component 1, q  
subcomponent binding protein

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A9B2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50719-1, C1QBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68084-2-PBS. Detection antibody: 68084-3-PBS. Standard:Ag19773. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50719-2, C1QBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68084-4-PBS. Detection antibody: 68084-3-PBS. Standard:Ag19773. Range: 0.195-25 ng/mL.