

For Research Use Only

C1QBP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68084-5-PBS



Basic Information

Catalog Number:

68084-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG19773

GenBank Accession Number:

BC013731

GeneID (NCBI):

708

UNIPROT ID:

Q07021

Full Name:

complement component 1, q
subcomponent binding protein

Calculated MW:

282 aa, 31 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E6A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50719-3, C1QBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68084-5-PBS. Detection antibody: 68084-6-PBS. Standard:Ag19773. Range: 0.391-100 ng/mL.