

For Research Use Only

FKBP10 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68085-2-PBS



Basic Information

Catalog Number: 68085-2-PBS	GenBank Accession Number: BC016467	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 60681	CloneNo.: 1F12A11
Source: Mouse	UNIPROT ID: Q96AY3	
Isotype: IgG1	Full Name: FK506 binding protein 10, 65 kDa	
Immunogen Catalog Number: AG31628	Calculated MW: 582 aa, 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50974-1, FKBP10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68085-2-PBS. Detection antibody: 68085-3-PBS. Standard:Ag31628. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50974-2, FKBP10 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68085-2-PBS. Detection antibody: 68085-4-PBS. Standard:Ag31628. Range: 0.195-100 ng/mL