

For Research Use Only

BNIP3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68091-2-PBS



Basic Information

Catalog Number: 68091-2-PBS	GenBank Accession Number: BC009342	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 664	CloneNo.: 1B8H6
Source: Mouse	UNIPROT ID: Q12983	
Isotype: IgG2b	Full Name: BCL2/adenovirus E1B 19kDa interacting protein 3	
Immunogen Catalog Number: AG27935	Calculated MW: 194 aa, 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50666-1, BNIP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68091-2-PBS. Detection antibody: 68091-3-PBS. Standard: Ag27935. Range: 0.098-100 ng/mL.