

For Research Use Only

MIST1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68092-3-PBS



Basic Information

Catalog Number:

68092-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18509

GenBank Accession Number:

BC113394

GeneID (NCBI):

168620

UNIPROT ID:

Q7RTS1

Full Name:

basic helix-loop-helix family,
member a15

Calculated MW:

189 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B1A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51171-1, BHLHA15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68092-2-PBS. Detection antibody: 68092-3-PBS. Standard:Ag18509. Range: 0.391-100 ng/mL.