

For Research Use Only

# SDCBP Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68096-3-PBS



## Basic Information

<b>Catalog Number:</b> 68096-3-PBS	<b>GenBank Accession Number:</b> BC113674	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6386	<b>CloneNo.:</b> 1G5G5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O00560	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> syndecan binding protein (syntenin)	
<b>Immunogen Catalog Number:</b> AG18078	<b>Calculated MW:</b> 298 aa, 32 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51226-1, SDCBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68096-2-PBS. Detection antibody: 68096-3-PBS. Standard:Ag18078. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP51226-2, SDCBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68096-3-PBS. Detection antibody: 68096-4-PBS. Standard:Ag18078. Range: 0.781-100 ng/mL.