

For Research Use Only

p115, USO1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68100-4-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 68100-4-PBS | GenBank Accession Number: BC032654 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 8615 | CloneNo.: 3D11B2 |
| Source: Mouse | UNIPROT ID: O60763 | |
| Isotype: IgG1 | Full Name: USO1 homolog, vesicle docking protein (yeast) | |
| Immunogen Catalog Number: AG5543 | Calculated MW: 962 aa, 108 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51177-2, p115, USO1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68100-4-PBS. Detection antibody: 68100-3-PBS. Standard: Ag5543. Range: 0.195-50 ng/mL.